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Niet-steriel product

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Producent

EN - CDS protocol -Cleaning, Disinfection and sterilization protocol

0 Intended use

OTNI 13 Reusable Surgical Instrument Set for OI Surgery is exclusively intended to be used as reusable surgical instruments during osseointegration surgery on femur and/or tibia. This IFU applies to the cleaning, disinfection and sterilization of the instruments. For surgical use of the instruments we hereby refer to the Surgical Technique document; separately included and also downloadable from our website and specified below.

1 Introduction, warnings and precautions

This document contains instructions to clean, disinfect, sterilise, inspect and maintain the OTNI 13 Reusable Surgical Instruments for Osseointegration; based on requirements of EN ISO 17664: Sterilisation of medical devices. The document also provides general rules for preparation before use of these instruments and instructs for inspection of instruments, including a timely decision of replacing an instrument. The maximum number of reprocessing the OTNI 13 Reusable Surgical Instruments for Osseointegration is limited to 500 cycles. The following series of steps to prepare instruments for use or replacement are given:

2 On usage location

Eliminate rough contamination by use of absorbing tissue within two hours post-operatively. Rins instruments with current water or place instruments in adlehyde-free disinfection tub (see attach 1).

3 Transport to location of cleaning and sterilisation

Prevent mechanical wearance by separating heavy and vulnerable instruments from one another during transport. Limit the time prior to the cleaning procedure to a maximum of 30 minutes after surgery.

4 Preparing for cleaning: disassembling of assembled parts

Disassemble OTNI 1307 Taper Installer and OTNI 1310 Corer.

5 Precleaning

Eliminate rough contamination, soak in cleaning detergent, brush, check moving parts, rinse in current water, inspect and place on absorbing tissue. The precleaning can be skipped if immediate manual cleaning and disinfection is executed. In case of instruments that are seriously contaminated, pre-cleaning by use of ultrasonic cleaning is preferred. Required method: A cleaning tub that is big enough to completely submerge/soak the instruments in the liquid and a freshmade cleaning detergent for manual cleaning (see attach 1). The temperature shall be conform the instructions of the manufacturer but not more than 50 °C. Do not use any metal brushes or steel wool. Use personal protection according to the preferences of the supplier of the cleaning detergent (at minimum you should wear overall, gloves and face/eye protection). The use of fresh made cleared or sterile water for rinsing is recommended.

6.1 Automatic Cleaning and Disinfection

It is preferred to clean and disinfect by the automatic thermo-disinfection procedure. This is better reproducable and therefore more trustworthy. Preferably make us of a Thermo-disinfector with CE conform ISO 15883.

Make use of a programm for thermal disinfection (AO-level >3000). A program for chemical disinfection is not recommended. Load the thermo-disinfector, connect hollow parts to the rinse -canals. Place the instruments such, that hollow parts are not placed horizontally and make sure that blind holes are facing downward. Instruments that consist of multiple parts must be placed in open position. Start the cycle of the thermodisinfecter (at minimum a period of 5 minutes on 90 °C). Make use of cleaning detergent that is developed to be used for thermal disinfection procedures and that meets the criteria and concentrations in attach 1. The final rinsing procedure is preferably performed exclusively with fresh made and cleared water. Dry the instrument by use of compressed air of medical quality and make use of cleaning /lintfree disposal tissues (if necessary you can add 2 hours of after drying on a clean location or by heating in an oven on maximum 110°C.

6.2 Manual Cleaning and Disinfection

Soak the instruments in ultrasonic tub (minimum 15 minutes, frequency 20–25 kHz. Use concentration and temperature according to the prescriptions of the manufacturer of the cleaning detergent. Make use of cleaning detergent that is meant to be used for manual cleaning (attach 1). Clean instruments by use of brushes (no metal brushes or steel wool) or cleaning threads. Make use of 50 cc injection syringes for a minimum of 3 times rinsing of canals, hollow parts and blind holes. Rinse instruments at minum 1 minute by use of spray / running water. Dry instruments by use of compressed air of medical quality and make use of cleaning / lint-free disposal tissues (if necessary you can add 2 hours of after drying on a clean location or by heating up in an oven on maximum 110 degrees Celsius.).

7 Inspection

Prior to sterilisation, inspect all instruments visually and check instruments for proper functioning. In case after the cleaning procedure some stick residu of contamination is still present on instruments, repeat the cleaning procedure. Check for contamination, corrosion, damage and wearance and check cutting surfaces of sharp /cutting instruments (corers,

burs and enlargers). It is recommended to have the OTNI 1310 Corer periodically sharpened; preferably prior to each surgical procedure. In case of damage or bluntness of instruments, contact the instrument department of the hospital /clinic or contact OTN Implants for repair or replacement.

8 Sterilization

Steam sterilization is recommended according to EN ISO 17665. The hospital will then need to validate their process of sterilization themselves. Cycle: saturated steam with fractional forced air-vacuum. Exposure time: 3 minutes. Temperature: 134–138° C. Drying time: minimum 30 minutes.

9 Drying

It is the responsibility of the hospital to determine and validate the suitable drying time for the machine that is used for the sterilization.

10 Packaging

If no instrument net is used during themachinal cleaning /disinfection and sterilization, it is recommended to make use of a double pouch. The pouch should meet the requirements of EN ISO 11607 and be suitable for steam sterilization. The pouch will need to offer sufficient protection of the instruments against mechanical wearance.

11 Storage prior to use

Check the environment, check the storage period. Store the medical devices after sterilization in the sterilization pouch on a dust free location. The storage periode is dependant from the packaging and the method of storage. The maximum storage period prior to use is to be determined by the hospital.

12 Manufacturer contact details

OTN Implants BV, Simon Stevinweg 48, 6827BT Arnhem, Netherlands
mail@otnimplants.nl, www.otnimplants.nl, telephone: +31 26 7370371.

13 References

- EN ISO 11607 Packaging of sterilized medical devices
- EN ISO 17665 Sterilization of products for health
- EN ISO 17664 Sterilization of medical devices. Information to be supplied by manufacturer for processing of resterilizable medical devices
- EN ISO 15883-5Wash-/disinfection machines

Attach 1

For cleaning and/or disinfection of medical instruments that have been produced of aluminum alloy, a pH neutral cleaning detergent is to be used. For automatic cleaning / disinfection machines, the use of a mild alcalic cleaning detergent is recommended (a pH of 10.9). It's recommended to exclusively make use of cleaning and disinfection detergents bearing CE marking. Make sure to comply to all instructions of the supplier of the cleaning detergent.

Recommended means to preclean and manually clean:

Cidezyme/Enzol (Johndson and Johnson).

Recommended means for manual disinfection:

Cidex OPA (Johnson and Johnson).

Recommended means for automatic cleaning/disinfection by use of a thermodisinfecter machine:

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